



Everlasting Beauty

Be successful with the
Permanent Make-up
system by LCN

Beauty Concepts International Ltd.
Danworth Farm, Hurstpierpoint
West Sussex BN6 9GL
Phone: 01273-837 302
Email: bci@beautyconcepts.co.uk
www.beautyconcepts.co.uk

Art.-Nr.: 635593-GB; 02/2010



Put your trust in a reliable partner

LCN has been a reliable partner in the field of nail and foot care.

By developing Permanent Make-up, we have created an additional beauty system with which you expand your service offers. Furthermore, permanent make-up provides another lucrative business opportunity.

Profit from satisfied customers

Permanent Make-up remains the shooting star of all cosmetic services. This is due to the enormous advantages of permanent make-up. A perfect look around the clock, beauty without effort, permanent good looks anytime, anywhere.

Let your institute profit from this enormous demand. Get your qualifications now in order to increase your success in your field of business. We offer the necessary training as well as colours and equipment, which ensure a safe and quality conscious treatment.



Success can be planned by combining safe and high quality products with competence.

Every woman wants to be beautiful, to be admired. This is the reason why many women seek professional advice from beauticians or a beauty salon of their choice. Your salon can fulfill this desire for permanent beauty: with the Permanent Make-up System by LCN.

Using tradition to achieve success

LCN stands for high quality and safe products as well as competence all that since 1914. In the course of its successful business development, our longstanding company has become one of the leading suppliers in the special fields of cosmetics and professional corrective systems in foot care.

LCN has also established itself in the innovative field Permanent Make-up. Today we are one of the market leaders worldwide for highly specialized high-tech methods for skin pigmentation.

The basis of LCN is inhouse research and production. Rigorous quality control is a must. Our employees cooperate with leading universities and specialists in the fields of research and development, consulting and of course service. The training standards are very high. These qualities sum up our company philosophy.



Let dreams come true

The desire for permanent/everlasting beauty is as old as mankind. Let a dream come true with the Permanent Make-up System by LCN.

Create permanent beauty

The principle of permanent make-up is as simple as it is effective: fine needles bring colour pigments into the skin. A few lines increase or discreetly correct Mother Nature's gifts permanently! Even small deficiencies can be covered up. The result: an attractive look, selfconfidence and success. Permanent Make-up – the revolution in the cosmetic business. Based on well founded scientific research and development, permanent make-up has been established as a perfectly safe, high-tech method for cosmetic treatment.



Pioneer equipment for excellent results

Reliable and perfect equipment provides the best basis for your success and inspires confidence in your customers. This is the guiding principle behind the continual research and innovative development of the precision instruments for permanent make-up by LCN.



Survey on the advantages of the Titanium unit

Hygiene module

- No contamination of applicator
- No autoclaving of parts
- No handling of the 'bare' needle
- Quick change of colours and change of needle combination
- Large choice of needle combinations
- Optimised needle properties
- Gamma sterilized, packaging in the clean room, Batch tracking

Applicator

- Extreme precision for finest lines
- Adjustable insertion depth
- Straight and slanted applicator available

Control unit

- Use of only high quality electronic parts
- Controlled by microprocessor
- Foot pedal
- Hygienic plastic foil keyboard

Persuasive technology – the Titanium unit

This Titanium unit with its innovative disposable hygiene module is the product developed along the lines of our guiding principle. Needle and nozzle are combined in a module.

Attach the hygiene module to the applicator. A membrane screen covers the back of this module. This protects the applicator against contaminated fluids and air.

The hygiene module allows for a quick exchange of needles or colours. When the hygiene module is removed from the applicator, the needle disappears immediately into the nozzle. The danger of handling contaminated needles and the imminent risk of infection (for example AIDS, hepatitis) are eliminated. The hygiene module is directly filled with the pigmenting colour.

Drawing superfine lines – maximum hygiene

The Titanium unit can be used for a wide variety of treatments, especially for drawing superfine lines. A microprocessor ensures the necessary regular insertion frequency. The Titanium unit can also be used for pigmenting larger areas (including the treatment of scars).

The control unit also has been designed with hygiene in mind. The unit is equipped with a plastic foil keyboard, which makes cleaning easy.

The Kryo-method – New technology for pigmentation

The latest research has produced a new type of unit by LCN; the Kryo-Liner. This high quality pigmenting unit has been equipped with a unique icing technique, which allows for an almost pain-free treatment using only a minimum of local anaesthetics. It reduces irritations to a minimum.

Cool and soft – the Kryo-Liner

The technical innovation of the Kryo-Liner consists of an integrated icing probe, which continuously indicates the current probe temperature. By correctly cooling the respective skin surface, the skin is optimally prepared for pigmentation.

Gentle treatment – precise pigmentation

The icing technique has several advantages. It almost eliminates the need for anaesthetics like creams or salves. Swelling and irritations of the treated skin areas are reduced due to the cooling. It prevents blood and lymph fluid from flowing. The skin absorbs the pigments.



Survey of the advantages of the Kryo-Liner

- Prevents micro-injuries of the skin
- Reduced swellings
- It causes vascular constriction, which prevents blood drops and lymph fluid to rise.
- Almost eliminates the need for local anaesthetics like creams or solutions
- Precise colouring through exact pigmentation
- Very good pigmenting results by reducing crusts and loss of pigments
- Clear lines due to less swellings and better preliminary drawings



Permanent Make-up Colours by LCN – How to enhance beauty

Tolerability, brilliance and the intensity of the colours have been a top priority for the development of Permanent Make-up Colours.

Dermatological tolerance

The assessment of the dermatological tolerance by a renowned institute of a German university comes to the following conclusions:

„that Permanent Make-up Colours regarding the different combinations of pigments in the colours can be considered safe.“

Biological safety

Internationally renowned scientific institutes for biological safety and dermatological tolerance have tested all ingredients.

Superior application properties

Permanent Make-up Colours are distinguished by the following properties:

- Very good adhesion to and application with the needle
- High brilliance of the colours due to the fine distribution of the pigments
- No skin reaction

There are no health risks

Evaluations done by the independent Swiss RCC Institute has been proven the products are non-toxic.

“Based on the respective data and calculations, a health risk while using the products assessed in this expertise can be excluded if they are applied correctly. The formula of the permanent make up colours ‘Permanent Make-up Colours’ do not contain components, which indicate a risk to one’s health, if applied correctly.

This expertise meets all requirements of the EC guideline 93/35/EWG on cosmetic products, section 1, paragraph 12 (subparagraph 7a (1) d).”

Ingredients

Permanent Make-up Colours consist of fine and pure pigments in a liquid. The colour of the pigments remains stable. The quality and the ingredients of the fluid and the pigments are described as follows:

- a) **Liquid:** Mixture of isopropyl alcohol, and physiological saline and glycerine. The quality and purity of the components meet the requirements of the German pharmacopoeia (DAB) as well as the US-pharmacopoeia (USP).
- b) **Pigments:** Blend of iron oxide pigments^{1,2}, titanium dioxide pigment^{1,2} and other special pigments³. The quality and purity of the pigments meet the standards and the European regulations for dyeing substances for drugs and cosmetic products.
- c) **Preservative:** Benzoic acid according the German regulations on cosmetic products³.

Permanent Make-up Special training. A stepping-stone to your success!

Lay the foundations for your success with a special training as a permanent make-up specialist! Profit from our experience and expertise in application, equipment, marketing and consulting.



Become a Permanent Make-up specialist

LCN training is not only aimed at beginners but also at advanced trainees. You will learn the latest techniques on permanent make-up and its application. We offer success-oriented methodical treatment. Experienced Permanent Make-up trainers hold the training seminars in authorized training centres. The training includes theoretical as well as practical training

Attract your customers with your certificate

After your successful exam you are awarded the Permanent Make-up certificate. This authorizes you to use this brand name for advertisement. You also have the opportunity to utilize practical training days for your advanced training.



For more information call our Customer Service team:

Make use of our expertise! Please don't hesitate to contact us for detailed information about beauty matters.
Call 01 27 3/83 73 02 or send us an
Email: bci@beautyconcepts.co.uk.